

# JAS NGA

Natural Gas Analyzer

jas Analyzer



# JAS NGA

Natural Gas Analyzer


**jas** Analyzer

**JAS NGA – a prime turnkey system for determination of natural gas composition and for calculation of BTU and gas properties.**

JAS NGA is based on the **Agilent's 8890 GC** that is configured for analysis of natural gas. With the combination of unique JAS isothermal valve box, valve automation technology, optimized GC column set and JAS Report Generator software, JAS NGA distinguishes itself from competitors as a prime turn-key system for natural gas analysis.

Included in the system, the JAS Report Generator is fully integrated into Agilent OpenLAB™ for BTU and gas property calculations based on ISO 6976.

**JAS NGA can be specially configured upon request to the following standard methods to ensure suitability for use in downstream processes.& methods:**



Method	Remarks
ASTM D1945	
ASTM D3588	NGA is your TURNKEY
GPA 2261	Available as 1- & 3-valve system
GPA 2286	Custom made upon request
GPA 2186	Custom made upon request
ISO 6974	Custom made upon request

## NGA is your TURNKEY SOLUTION for natural gas analysis

- Turnkey solution — guaranteed application as custom-made solution to analyze one of followings: Natural Gas, Liquefied Petroleum Gas or Natural Gas Liquids (NGL)
- Automatic calculation of gas properties (e.g. BTU, specific gravity) and report generation at the end of an analysis (fully integrated into the Agilent OpenLAB)
- Operational procedures are fully documented – the analyzer not only incorporates proven GC hardware and software, but arrives also with the pre-loaded analysis method(s), documentation and a calculation model specific to the application
- Flexibility to analyze natural gas, liquefied petroleum gas or natural gas liquids (NGL)
- Rugged and easy to maintain
- Configuration alternatively possible with Agilent 8860 GC
- Configuration according to GPA 2261 and ASTM D 1945 possible; GPA 2186, GPA 2286, ISO 6974 upon request

## FIELD of APPLICATIONS

- Natural gas plants
- Pipeline monitoring stations
- Gas-fired power plants
- Gas Companies

# JAS NGA

Natural Gas Analyzer

jas Analyzer

## JAS REPORT GENERATOR, A DEDICATED AND UNIQUE SOFTWARE

This Software has the capability to automatically calculate physical properties such as BTU and specific gravity selected by the user after each analysis. The S/W uses the result table generated in OpenLAB™ to calculate needed properties according to ISO 6976. A standard database for gases with flexibility for users to modify for special needed is included in the software. Additional conditions such as temperature and pressure can be edited as well.

## SOFTWARE FEATURES

- Expandable Compound Table
- Calculations according to ISO 6976
- Flexible Report Generation, allows customization

The screenshot displays the JAS BTUcalc software interface. It features a table of compounds with columns for CalNr, Compound, MolWeight (g/mol), Summation Factor, and Gross Calorific Value (MJ/m³). The table lists 17 compounds, including Hydrogen, Carbon Monoxide, Carbon Dioxide, Oxygen, Nitrogen, Hydrogensulfide, Ethylene, Methane, Acetylene, Propylene, Propadiene, Methylacetylene, and t-2-Butene. To the right of the table, there are input fields for 'Metering constants' (Temperature (K), Pressure (Pa), Mass of Air (KG/K.mol), Compression Factor) and 'Calculations' (Compressibility (Zmix), Superior Calorific Value on volumetric, mass, and molar basis, Relative Density of Real Gas, Density of Real Gas (Kg/m³), Wobbe Index).

CalNr	Compound	MolWeight g/mol	Summation Factor	Gross Calorific Value MJ/m³
01	Hydrogen	2,0159	-0,0048	12,1020
04	Carbon Monoxide	28,0104	0,0224	11,9600
05	Carbon Dioxide	44,0098	0,0748	0,0000
06	Oxygen	31,9998	0,0283	0,0000
07	Nitrogen	28,0135	0,0173	0,0000
08	Hydrogensulfide	34,0819	0,1000	23,7800
09	Ethylene	28,0540	0,0800	59,7200
10	Methane	16,0428	0,0447	37,7060
11	Acetylene	26,0380	0,0837	55,0400
12	Propylene	42,0810	0,1265	87,1000
13	Propadiene	40,0648	0,1304	82,2100
14	Methylacetylene	40,0648	0,1304	82,2100
17	t-2-Butene	56,1080	0,1789	114,5400

**Metering constants**

Temperature (K) = 288,15  
Pressure (Pa) = 101,325  
Mass of Air (KG/K.mol) = 28,9626  
Compression Factor = 0,99958

**Calculations**

Compressibility (Zmix) = 0,9766318895  
Superior Calorific Value on volumetric basis = 102,7149921285 MJ/m³ (H)  
Superior Calorific Value on mass basis = 49,4898561052 MJ/Kg  
Superior Calorific Value on molar basis = 2.428,6882770891 KJ/mol  
Relative Density of Real Gas = 1,7336526770  
Density of Real Gas (Kg/m³) = 2,1253274592 Kg/m³  
Wobbe Index = 78,0104512252 MJ/m³

# JAS NGA

Natural Gas Analyzer

## INSTRUMENT CONFIGURATION according to GPA 2261

Hardware: Agilent GC 8890, TCD, FID, Valves, PC Bundle, Capillary and Packed Columns, Unis-S/SL, OpenLAB


Software: Special Calculation & Reporting Software functions, fully integrated into Agilent's OpenLAB

## Natural Gas COMPONENTS

The below table shows the natural gas components and their concentration range in %.

NGA-components	Formula	Concentration range in %
Oxygen	O <sub>2</sub>	0.01 ~ 20
Nitrogen	N <sub>2</sub>	0.01 ~ 100
Carbon Dioxide	CO <sub>2</sub>	0.01 ~ 20
Methane	CH <sub>4</sub>	0.01 ~ 100
Ethane	C <sub>2</sub> H <sub>6</sub>	0.01 ~ 15
Hydrogen sulphite	H <sub>2</sub> S	3.00 ~ 100
Propane	C <sub>3</sub> H <sub>8</sub>	0.01 ~ 100
iso-Butane	C <sub>4</sub> H <sub>10</sub>	0.01 ~ 10
n-Butane	C <sub>4</sub> H <sub>10</sub>	0.01 ~ 10
iso-Pentane	C <sub>5</sub> H <sub>12</sub>	0.01 ~ 2
n-Pentane	C <sub>5</sub> H <sub>12</sub>	0.01 ~ 2
Hexane and heavier components	C <sub>6</sub> H <sub>14</sub>	0.01 ~ 2

## Chromatogram



The following chromatogram shows the components C<sub>3</sub>-C<sub>5</sub>, CO<sub>2</sub>, C<sub>2</sub>H<sub>6</sub> as well as O<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub> and C<sub>6+</sub>. Back-flushed; based on 1-valve GPA 2261 analysis

# JAS NGA


Natural Gas Analyzer

**jas** Analyzer

## OPTION GICU – GAS INJECTION CONTROL UNIT

The GICU System is a standalone modular system and fully controlled by dedicated software that is fully integrated in Agilent OpenLAB™. The GICU allows gases to fill the GC sample loop with user specified pressures. By varying the pressure, the sample concentration can be controlled:

- As the pressure in the sample loop can be controlled with the GICU system, a multi-level calibration is possible with only using one standard gas.
- GICU allows sub ambient pressure injection of high concentration samples to avoid detector saturation and ensure samples are analyzed within the dynamic range.
- GICU allows constant pressure in the sampling loop and greatly increases the repeatability of gas sample injection.



## GICU BENEFITS:

- Very easy to calibrate and run gas samples under the same or similar conditions.
- GICU can handle different gas samples with different pressures
- Possibility to adjust the sample loop pressure from 50 mbar to 1800 mbar (abs) / 0.725 psi to 26 psi (abs.)
- GICU is using an automated sampling system with a pre-programmed procedure for flushing and cleaning of the gas path of the system: ideal when handling hazardous samples.
- The GICU System is modular and stand alone, fully controlled by dedicated software: GICU starts the GC run after achieving the set pressure or when the GICU method is complete.
- Setup of the program is easy and user-friendly and can be incorporated into specific selected GC methods



## About Joint Analytical Systems

Since 1995 JAS has been a Premier Solution Partner and Value Added Reseller of Agilent Technologies. We are an innovative-driven organization that offers customized solutions for GC, GCxGC,  $\mu$ GC, GC-AED, GC-MS, GC-QQQ, LC, LC-MS, LC-QQQ and Q-TOF LC-MS applications.

## JAS serves key industries such as

- Chemical
- Petrochemical/HPI
- Environmental
- Food & Flavor
- Forensic

## JAS Products for GC

- AED – Atomic Emission Detector
- Automatic Gas Samplers
- Customized Valving Systems
- EzPrep – Preparative Fraction Collector
- GICU – Gas Injection Control Unit
- UNIS Inlet Systems
- Olfactometer
- CryoTrap

## Joint Analytical Systems GmbH

Carl-Zeiss-Straße 49  
47445 Moers  
Germany

Phone: +49 2841 9871 100  
Fax: +49 2841 9871 222  
e-Mail: [info@jas.de](mailto:info@jas.de)  
Internet: [www.jas.de](http://www.jas.de)

Copyright © 2025  
Joint Analytical Systems GmbH  
All rights reserved